

Product datasheet for **TS423188P5**

PPAP2A (PLPP1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PLPP1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223188
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PPAP2A (PLPP1) Mouse Monoclonal Antibody [Clone ID: OT11H4] (TA506857)
ACCN:	NM_176895 , NP_795714
Synonyms:	LLP1a; LPP1; PAP-2a; PAP2; PPAP2A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_795714
Locus ID:	8611
Cytogenetics:	5q11.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Ether lipid metabolism, Fc gamma R-mediated phagocytosis, Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways, Sphingolipid metabolism



[View online »](#)