

Product datasheet for **TS423222**

beta Defensin 3 (DEFB103B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DEFB103B in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223222
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_018661 , NP_061131
Synonyms:	BD-3; DEFB-3; DEFB3; DEFB103; HBD-3; HBD3; HBP-3; HBP3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_061131
Locus ID:	55894
Cytogenetics:	8p23.1
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)