

Product datasheet for **TS423314P5**

Macrophage Scavenger Receptor I (MSR1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MSR1, transcript variant SR-AIII, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223314
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Macrophage Scavenger Receptor I (MSR1) Mouse Monoclonal Antibody [Clone ID: OT19E5] (TA802777)
ACCN:	NM_138716 , NP_619730
Synonyms:	CD204; phSR1; phSR2; SCARA1; SR-A; SR-AI; SR-AII; SR-AIII; SRA
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_619730
Locus ID:	4481
Cytogenetics:	8p22
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)