

## Product datasheet for **TS423438P5**

### RNF14 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RNF14, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223438
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RNF14 Mouse Monoclonal Antibody [Clone ID: OTI1C9] (TA811050)
ACCN:	<a href="#">NM_183398</a> , <a href="#">NP_899645</a>
Synonyms:	ARA54; HFB30; HRIHFB2038; TRIAD2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_899645</a>
Locus ID:	9604
Cytogenetics:	5q31.3
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)