

Product datasheet for **TS423486P5**

Angiogenin (ANG) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANG, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223486
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001097577</u> , <u>NP_001091046</u>
Synonyms:	ALS9; HEL168; RAA1; RNASE4; RNASE5
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001091046</u>
Locus ID:	283
Cytogenetics:	14q11.2
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)