

Product datasheet for **TS423510P5**

Eph receptor A6 (EPHA6) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EPHA6, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223510
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Eph receptor A6 (EPHA6) Mouse Monoclonal Antibody [Clone ID: OT15B9] (TA809373)
ACCN:	NM_001080448 , NP_001073917
Synonyms:	EHK-2; EHK2; EK12; EPA6; HEK12; PRO57066
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001073917
Locus ID:	285220
Cytogenetics:	3q11.2
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Axon guidance



[View online »](#)