

## Product datasheet for **TS423537**

### Exosome Component 9 (EXOSC9) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EXOSC9, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223537
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001034194</a></u> , <u><a href="#">NP_001029366</a></u>
Synonyms:	p5; p6; PCH1D; PM/ScI-75; PMSCL1; RRP45; Rrp45p
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001029366</a></u>
Locus ID:	5393
Cytogenetics:	4q27
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	RNA degradation



[View online »](#)