

## **Product datasheet for TS423561P5**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **MAGEA2 CytoSection**

**Product data:** 

**Product Type:** CytoSections

**Description:** Transient overexpression of MAGEA2, transcript variant 1, in HEK293T cells, FFPE control for

IHC, ICC and ISH staining, 25 slides per pack

Species: Human
Expression Host: HEK293T

**Expression cDNA Clone** 

or AA Sequence:

TrueORF Clone RC223561

Tag: C-MYC/DDK

**Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

ACCN: <u>NM 005361, NP 005352</u>

Synonyms: CT1.2; MAGE2; MAGEA2A

Storage: Room Temperature

**Stability:** Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

**Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

**Note:** This product is for research use only and is not approved for use in humans or in clinical


diagnosis.

RefSeq: NP 005352


Locus ID: 4101 Cytogenetics: Xq28



## **Product images:**



Immunohistochemistry staining on MAGEA4 variant 4 overexpressed cytosection TS423561P5 with mouse anti DDK clone OTI11C3 C/N [TA180144] at 1:600 dilution 1hr room temperature. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D12-18]). Positive stain shown with the brown chromogen present. mammalian cells transfected with cDNA clone [RC223561], 5 micron sections, 40x magnification.



Immunohistochemistry staining on MAGEA2 overexpressed cytosection TS423561P5 with mouse anti DDK clone OTI11C3 C/N [TA180144] at 1:600 dilution 1hr room temperature. Antibody staining was achieved with HIER Citrate pH6, Polink1 with DAB chromogen detection kit ([D12-18]). Positive stain shown with the brown chromogen present. mammalian cells transfected with cDNA clone [RC223561], 5 micron sections, 40x magnification.