

Product datasheet for **TS423612P5**

Haptoglobin (HP) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HP, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223612
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Haptoglobin (HP) Mouse Monoclonal Antibody [Clone ID: OTI4H5] (TA501699)
ACCN:	NM_005143 , NP_005134
Synonyms:	BP; HP2ALPHA2; HPA1S
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005134
Locus ID:	3240
Cytogenetics:	16q22.2
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane



[View online »](#)