

## Product datasheet for **TS423691P5**

### TAF4 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TAF4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223691
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_003185</a> , <a href="#">NP_003176</a>
Synonyms:	TAF(II)130; TAF(II)135; TAF2C; TAF2C1; TAF4A; TAFII-130; TAFII-135; TAFII130; TAFII135
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_003176</a>
Locus ID:	6874
Cytogenetics:	20q13.33
Protein Pathways:	Basal transcription factors, Huntington's disease



[View online »](#)