

## Product datasheet for **TS423979**

### HLA-DQA2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HLA-DQA2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223979
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HLA-DQA2 Mouse Monoclonal Antibody [Clone ID: OTI4C9] (TA812936)
ACCN:	<a href="#">NM_020056</a> , <a href="#">NP_064440</a>
Synonyms:	DC-alpha; DX-ALPHA; HLA-DCA; HLA-DXA; HLA-DQA2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_064440</a>
Locus ID:	3118
Cytogenetics:	6p21.32
Protein Families:	Transmembrane



[View online »](#)

**Protein Pathways:**

Allograft rejection, Antigen processing and presentation, Asthma, Autoimmune thyroid disease, Cell adhesion molecules (CAMs), Graft-versus-host disease, Systemic lupus erythematosus, Type I diabetes mellitus, Viral myocarditis