

Product datasheet for **TS424067**

PIG3 (TP53I3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TP53I3, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224067
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PIG3 (TP53I3) Mouse Monoclonal Antibody [Clone ID: OTI3B11] (TA503656)
ACCN:	<u>NM_147184</u> , <u>NP_671713</u>
Synonyms:	PIG3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_671713</u>
Locus ID:	9540
Cytogenetics:	2p23.3
Protein Families:	Druggable Genome
Protein Pathways:	p53 signaling pathway



[View online »](#)