

## Product datasheet for **TS424189**


### Histone H1.5 (HIST1H1B) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of H1-5 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224189
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_005322</a> , <a href="#">NP_005313</a>
Synonyms:	H1; H1.5; H1B; H1F5; H1s-3; HIST1H1B
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_005313</a>
Locus ID:	3009
Cytogenetics:	6p22.1



[View online »](#)

**Product images:**

Immunohistology staining on HIST1H1B overexpressed cytosection TS424189 with mouse anti DDK clone OT111C3 C/N [TA180144] at 1:400 dilution 4C O/N. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone [RC224189], 5 micron sections, 40x magnification