

## Product datasheet for **TS424398P5**

### ErbB 3 (ERBB3) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ERBB3, transcript variant s, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224398
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ErbB 3 (ERBB3) Mouse Monoclonal Antibody [Clone ID: OT13E10] (TA803694)
ACCN:	<a href="#">NM_001005915</a> , <a href="#">NP_001005915</a>
Synonyms:	c-erbB-3; c-erbB3; ErbB-3; erbB3-S; FERLK; HER3; LCCS2; MDA-BF-1; p45-sErbB3; p85-sErbB3; p180-ErbB3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001005915</a>
Locus ID:	2065
Cytogenetics:	12q13.2
Protein Families:	Adult stem cells, Druggable Genome, Protein Kinase, Secreted Protein, Stem cell - Pluripotency, Transmembrane



[View online »](#)

**Protein Pathways:** Calcium signaling pathway, Endocytosis, ErbB signaling pathway