

Product datasheet for **TS424569**

Lamin B Receptor (LBR) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LBR, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224569
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_194442 , NP_919424
Synonyms:	C14SR; DHCR14B; LMN2R; PHA; PHASK; TDRD18
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_919424
Locus ID:	3930
Cytogenetics:	1q42.12
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)