

## Product datasheet for **TS424767**

### Neurogenin3 (NEUROG3) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NEUROG3 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224767
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Neurogenin3 (NEUROG3) Mouse Monoclonal Antibody [Clone ID: OTI3B5] (TA500346)
ACCN:	<a href="#">NM_020999</a> , <a href="#">NP_066279</a>
Synonyms:	Atoh5; bHLHa7; Math4B; NGN-3; ngn3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_066279</a>
Locus ID:	50674
Cytogenetics:	10q22.1
Protein Families:	ES Cell Differentiation/IPS
Protein Pathways:	Maturity onset diabetes of the young



[View online »](#)