

Product datasheet for **TS424784P5**

IL36 gamma (IL36G) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IL36G in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224784
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	IL36 gamma (IL36G) Mouse Monoclonal Antibody [Clone ID: OTI2F4] (TA505994)
ACCN:	NM_019618 , NP_062564
Synonyms:	IL-1F9; IL-1H1; IL-1RP2; IL1E; IL1F9; IL1H1; IL1RP2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_062564
Locus ID:	56300
Cytogenetics:	2q14.1
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)