

Product datasheet for **TS424830**

BNIP1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of BNIP1, transcript variant BNIP1-b, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224830
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	BNIP1 Mouse Monoclonal Antibody [Clone ID: OTI2G8] (TA802217)
ACCN:	<u>NM_013979</u> , <u>NP_053582</u>
Synonyms:	NIP1; SEC20; TRG-8
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_053582</u>
Locus ID:	662
Cytogenetics:	5q35.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport



[View online »](#)