

Product datasheet for **TS424904**

SynCAM (CADM1) CytoSection

Product data:

| | |
|---------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CytoSections |
| Description: | Transient overexpression of CADM1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC224904 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | NM_001098517 , NP_001091987 |
| Synonyms: | BL2; IGSF4; IGSF4A; Necl-2; NECL2; RA175; sgiGSF; ST17; sTSLC-1; SYNCAM; synCAM1; TSLC1 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_001091987 |
| Locus ID: | 23705 |
| Cytogenetics: | 11q23.3 |
| Protein Families: | Transmembrane |
| Protein Pathways: | Cell adhesion molecules (CAMs) |



[View online »](#)