

Product datasheet for **TS424928**

Nogo B receptor (NUS1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NUS1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224928
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Nogo B receptor (NUS1) Mouse Monoclonal Antibody [Clone ID: OTI6C11] (TA803142)
ACCN:	NM_138459 , NP_612468
Synonyms:	C6orf68; CDG1AA; MGC:7199; MRD55; NgBR; TANGO14
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_612468
Locus ID:	116150
Cytogenetics:	6q22.1



[View online »](#)