

## Product datasheet for **TS425085**

### CD59 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CD59, transcript variant 7, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225085
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD59 Mouse Monoclonal Antibody [Clone ID: OTI2E11] (TA507349)
ACCN:	<a href="#">NM_001127226</a> , <a href="#">NP_001120698</a>
Synonyms:	1F5; 16.3A5; EJ16; EJ30; EL32; G344; HRF-20; HRF20; MAC-IP; MACIF; MEM43; MIC11; MIN1; MIN2; MIN3; MIRL; MSK21; p18-20
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001120698</a>
Locus ID:	966
Cytogenetics:	11p13
Protein Families:	Druggable Genome
Protein Pathways:	Complement and coagulation cascades, Hematopoietic cell lineage



[View online »](#)