

Product datasheet for **TS425332P5**

CD16 (FCGR3A) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FCGR3A, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225332
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD16 (FCGR3A) Mouse Monoclonal Antibody [Clone ID: OT12A2] (TA813632)
ACCN:	NM_001127593 , NP_001121065
Synonyms:	CD16; CD16A; FCG3; FCGR3; FCGR3III; FCR-10; FCR3III; FCR3IIIA; IGFR3; IMD20
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001121065
Locus ID:	2214
Cytogenetics:	1q23.3
Protein Families:	ES Cell Differentiation/IPS, Secreted Protein, Transmembrane
Protein Pathways:	Fc gamma R-mediated phagocytosis, Natural killer cell mediated cytotoxicity, Systemic lupus erythematosus



[View online »](#)