

## Product datasheet for **TS425423P5**

### Claudin 5 (CLDN5) CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of CLDN5, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack   |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC225423  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| Target Detection Antibodies:          | Claudin 5 (CLDN5) Mouse Monoclonal Antibody [Clone ID: OTI1G4] (TA500843)   |
| ACCN:                                 | <a href="#">NM_001130861</a> , <a href="#">NP_001124333</a>   |
| Synonyms:                             | AWAL; BEC1; CPETRL1; TMDVCF; TMVCF  |
| Storage:                              | Room Temperature  |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <a href="#">NP_001124333</a>  |
| Locus ID:                             | 7122  |
| Cytogenetics:                         | 22q11.21  |
| Protein Families:                     | Transmembrane   |
| Protein Pathways:                     | Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction  |



[View online »](#)