

## Product datasheet for **TS425607**

### p63 (TP63) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TP63, transcript variant 6, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225607
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	p63 (TP63) Mouse Monoclonal Antibody [Clone ID: OTI4E7] (TA802078)
ACCN:	<u><a href="#">NM_001114982</a></u> , <u><a href="#">NP_001108454</a></u>
Synonyms:	AIS; B(p51A); B(p51B); EEC3; KET; LMS; NBP; OFC8; p40; p51; p53CP; p63; p73H; p73L; RHS; SHFM4; TP53CP; TP53L; TP73L
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001108454</a></u>
Locus ID:	8626
Cytogenetics:	3q28
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)