

Product datasheet for **TS425693P5**

ACADM CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACADM, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225693
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ACADM Mouse Monoclonal Antibody [Clone ID: OTI2G7] (TA811761)
ACCN:	<u>NM_001127328</u> , <u>NP_001120800</u>
Synonyms:	ACAD1; MCAD; MCADH
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001120800</u>
Locus ID:	34
Cytogenetics:	1p31.1
Protein Families:	Druggable Genome
Protein Pathways:	beta-Alanine metabolism, Fatty acid metabolism, Metabolic pathways, PPAR signaling pathway, Propanoate metabolism, Valine, leucine and isoleucine degradation



[View online »](#)