

Product datasheet for **TS425700**

Cbl c (CBLC) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CBLC, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225700
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Cbl c (CBLC) Mouse Monoclonal Antibody [Clone ID: OTI3D9] (TA505052)
ACCN:	NM_001130852 , NP_001124324
Synonyms:	CBL-3; CBL-SL; RNF57
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001124324
Locus ID:	23624
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways:

Chronic myeloid leukemia, Endocytosis, ErbB signaling pathway, Insulin signaling pathway, Jak-STAT signaling pathway, Pathways in cancer, T cell receptor signaling pathway, Ubiquitin mediated proteolysis