

Product datasheet for **TS425988P5**

MSF (SEPT9) CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of SEPTIN9, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC225988 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | MSF (SEPT9) Mouse Monoclonal Antibody [Clone ID: OTI13B1] (TA812821) |
| ACCN: | NM_001113491 , NP_001106963 |
| Synonyms: | AF17q25; MSF; MSF1; NAPB; PNUTL4; SEPT9; SeptD1; SINT1 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_001106963 |
| Locus ID: | 10801 |
| Cytogenetics: | 17q25.3 |
| Protein Families: | Druggable Genome |



[View online »](#)