

## Product datasheet for **TS426017P5**

### MMP2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MMP2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC226017
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MMP2 Mouse Monoclonal Antibody [Clone ID: OTI4A11] (TA806846)
ACCN:	<a href="#">NM_001127891</a> , <a href="#">NP_001121363</a>
Synonyms:	CLG4; CLG4A; MMP-2; MMP-II; MONA; TBE-1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001121363</a>
Locus ID:	4313
Cytogenetics:	16q12.2
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Bladder cancer, GnRH signaling pathway, Leukocyte transendothelial migration, Pathways in cancer



[View online »](#)