

Product datasheet for **TS426147**

CD71 (TFRC) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TFRC, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC226147
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD71 (TFRC) Mouse Monoclonal Antibody [Clone ID: OT16H9] (TA803640)
ACCN:	NM_001128148 , NP_001121620
Synonyms:	CD71; IMD46; p90; T9; TFR; TFR1; TR; TRFR
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001121620
Locus ID:	7037
Cytogenetics:	3q29
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Protease, Secreted Protein, Transmembrane
Protein Pathways:	Endocytosis, Hematopoietic cell lineage



[View online »](#)