

## Product datasheet for **TS426279P5**

### **PAG3 (ASAP2) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of ASAP2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC226279
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>ACCN:</b>	<u><a href="#">NM_001135191</a></u> , <u><a href="#">NP_001128663</a></u>
<b>Synonyms:</b>	AMAP2; CENTB3; DDEF2; PAG3; PAP; Pap-alpha; SHAG1
<b>Storage:</b>	Room Temperature
<b>Stability:</b>	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u><a href="#">NP_001128663</a></u>
<b>Locus ID:</b>	8853
<b>Cytogenetics:</b>	2p24
<b>Protein Pathways:</b>	Endocytosis, Fc gamma R-mediated phagocytosis



[View online »](#)