

Product datasheet for **TS426587**

NDUFB11 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NDUFB11, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC226587
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NDUFB11 Mouse Monoclonal Antibody [Clone ID: OT19C5] (TA807599)
ACCN:	NM_001135998 , NP_001129470
Synonyms:	CI-ESSS; ESSS; MC1DN30; Np15; NP17.3; P17.3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001129470
Locus ID:	54539
Cytogenetics:	Xp11.3
Protein Families:	Transmembrane



[View online »](#)