

## Product datasheet for **TS426680**

### NFYC CytoSection

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | CytoSections   |
| Description:                          | Transient overexpression of NFYC, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack  |
| Species:                              | Human  |
| Expression Host:                      | HEK293T  |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC226680   |
| Tag:                                  | C-MYC/DDK  |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)   |
| ACCN:                                 | <u>NM_001142589</u> , <u>NP_001136061</u>  |
| Synonyms:                             | CBF-C; CBFC; H1TF2A; HAP5; HSM; NF-YC  |
| Storage:                              | Room Temperature   |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.  |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.  |
| RefSeq:                               | <u>NP_001136061</u>  |
| Locus ID:                             | 4802   |
| Cytogenetics:                         | 1p34.2   |
| Protein Families:                     | Transcription Factors  |
| Protein Pathways:                     | Antigen processing and presentation  |



[View online »](#)