

Product datasheet for **TS426697**

MYH14 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MYH14, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC226697
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001145809</u> , <u>NP_001139281</u>
Synonyms:	DFNA4; DFNA4A; FP17425; MHC16; MYH17; myosin; NMHC-II-C; NMHC II-C; PNMHH
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001139281</u>
Locus ID:	79784
Cytogenetics:	19q13.33
Protein Pathways:	Regulation of actin cytoskeleton, Tight junction, Viral myocarditis



[View online »](#)