

## **Product datasheet for TS427317P5**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **FGF1 CytoSection**

**Product data:** 

**Product Type:** CytoSections

Description: Transient overexpression of FGF1, transcript variant 4, in HEK293T cells, FFPE control for IHC,

ICC and ISH staining, 25 slides per pack

Species: Human **Expression Host:** HEK293T

**Expression cDNA Clone** 

or AA Sequence:

TrueORF Clone RC227317

C-MYC/DDK Tag:

**Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

**Target Detection** 

Antibodies:

FGF1 Mouse Monoclonal Antibody [Clone ID: OTI4D2] (TA500064)

ACCN: NM 001144892, NP 001138364

AFGF; ECGF; ECGF-beta; ECGFA; ECGFB; FGF-1; FGF-alpha; FGFA; GLIO703; HBGF-1; HBGF1 Synonyms:

Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

This product is for research use only and is not approved for use in humans or in clinical Note:

diagnosis.

RefSeq: NP 001138364

Locus ID: 2246 Cytogenetics: 5q31.3

**Protein Families:** Druggable Genome, Secreted Protein

**Protein Pathways:** MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton

