

## Product datasheet for **TS427358P5**

### Mortality Factor 4 like 2 (MORF4L2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MORF4L2, transcript variant 15, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227358
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001142431</u> , <u>NP_001135903</u>
Synonyms:	MORFL2; MRGX
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001135903</u>
Locus ID:	9643
Cytogenetics:	Xq22.2
Protein Families:	Transcription Factors



[View online »](#)