

Product datasheet for **TS427644P5**

SnoN (SKIL) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SKIL, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227644
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SnoN (SKIL) Mouse Monoclonal Antibody [Clone ID: OTI2F6] (TA800171)
ACCN:	<u>NM_001145097</u> , <u>NP_001138569</u>
Synonyms:	SNO; SnoA; Snol; SnoN
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001138569</u>
Locus ID:	6498
Cytogenetics:	3q26.2
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)