

Product datasheet for **TS427787P5**

EYS CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EYS, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227787
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001142801</u> , <u>NP_001136273</u>
Synonyms:	bA74E24.1; bA166P24.2; bA307F22.3; C6orf178; C6orf179; C6orf180; dj22i17.2; dj303f19.1; dj1018A4.2; EGFL10; EGFL11; RP25; SPAM
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001136273</u>
Locus ID:	346007
Cytogenetics:	6q12
Protein Families:	Druggable Genome



[View online »](#)