

## Product datasheet for **TS427813**

### CD299 (CLEC4M) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CLEC4M, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227813
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD299 (CLEC4M) Mouse Monoclonal Antibody [Clone ID: OTI7D12] (TA810067)
ACCN:	<a href="#">NM_001144911</a> , <a href="#">NP_001138383</a>
Synonyms:	CD209L; CD299; DC-SIGN2; DC-SIGNR; DCSIGNR; HP10347; L-SIGN; LSIGN
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001138383</a>
Locus ID:	10332
Cytogenetics:	19p13.2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)