

Product datasheet for **TS427975P5**

Ki67 (MKI67) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MKI67, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227975
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Ki67 (MKI67) Rat Monoclonal Antibody [Clone ID: OTI5D7] (TA801577)
ACCN:	NM_001145966 , NP_001139438
Synonyms:	KIA; MIB-; MIB-1; PPP1R105
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001139438
Locus ID:	4288
Cytogenetics:	10q26.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS



[View online »](#)