

Product datasheet for **TS428276P5**

Aminomethyltransferase (AMT) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AMT, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC228276
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Aminomethyltransferase (AMT) Mouse Monoclonal Antibody [Clone ID: OT16F4] (TA807784)
ACCN:	NM_001164711 , NP_001158183
Synonyms:	GCE; GCST; GCVT; NKH
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001158183
Locus ID:	275
Cytogenetics:	3p21.31
Protein Pathways:	Glycine, serine and threonine metabolism, Metabolic pathways, Nitrogen metabolism, One carbon pool by folate



[View online »](#)