

## Product datasheet for **TS428689P5**

### PEG3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PEG3, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC228689
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PEG3 Mouse Monoclonal Antibody [Clone ID: OTI1E10] (TA801960)
ACCN:	<a href="#">NM_001146186</a> , <a href="#">NP_001139658</a>
Synonyms:	PW1; ZKSCAN22; ZNF904; ZSCAN24
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001139658</a>
Locus ID:	5178
Cytogenetics:	19q13.43
Protein Families:	Transcription Factors



[View online »](#)