

## Product datasheet for **TS429072**

### **KDM3A / JHDM2A (KDM3A) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of KDM3A, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC229072
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	KDM3A / JHDM2A (KDM3A) Mouse Monoclonal Antibody [Clone ID: OT15A7] (TA812138)
<b>ACCN:</b>	<a href="#">NM_001146688</a> , <a href="#">NP_001140160</a>
<b>Synonyms:</b>	JHDM2A; JHMD2A; JMJD1; JMJD1A; TSGA
<b>Storage:</b>	Room Temperature
<b>Stability:</b>	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<a href="#">NP_001140160</a>
<b>Locus ID:</b>	55818
<b>Cytogenetics:</b>	2p11.2



[View online »](#)