

Product datasheet for **TS429677P5**

Galectin 3 (LGALS3) CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of LGALS3, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC229677 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | Galectin 3 (LGALS3) Mouse Monoclonal Antibody [Clone ID: OTI1G8] (TA506399) |
| ACCN: | <u>NM_001177388</u> , <u>NP_001170859</u> |
| Synonyms: | CBP35; GAL3; GALBP; GALIG; L31; LGALS2; MAC2 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_001170859</u> |
| Locus ID: | 3958 |
| Cytogenetics: | 14q22.3 |
| Protein Families: | Secreted Protein |



[View online »](#)