

Product datasheet for **TS429809P5**

CCRK (CDK20) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CDK20, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229809
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CCRK (CDK20) Mouse Monoclonal Antibody [Clone ID: OT11B3] (TA809548)
ACCN:	NM_001170639 , NP_001164110
Synonyms:	CCRK; CDCH; P42; PNQALRE
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001164110
Locus ID:	23552
Cytogenetics:	9q22.1
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)