

## Product datasheet for **TS429912**

### PPAR delta (PPARD) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PPARD, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229912
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PPAR delta (PPARD) Mouse Monoclonal Antibody [Clone ID: OTI4D10] (TA811238)
ACCN:	<a href="#">NM_001171820</a> , <a href="#">NP_001165291</a>
Synonyms:	FAAR; NR1C2; NUC1; NUCI; NUCII; PPARB
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001165291</a>
Locus ID:	5467
Cytogenetics:	6p21.31
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	Acute myeloid leukemia, Pathways in cancer, PPAR signaling pathway, Wnt signaling pathway



[View online »](#)