

Product datasheet for **TS430010P5**

Cell adhesion molecule 2 (CADM2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CADM2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230010
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001167675 , NP_001161147
Synonyms:	IGSF4D; Necl-3; NECL3; SynCAM 2; synCAM2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001161147
Locus ID:	253559
Cytogenetics:	3p12.1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)