

Product datasheet for **TS430106P5**

CP2c (TFCP2) CytoSection

Product data:

| | |
|---------------------------------------|--|
| Product Type: | CytoSections |
| Description: | Transient overexpression of TFCP2, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC230106 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | <u>NM_001173453</u> , <u>NP_001166924</u> |
| Synonyms: | LBP1C; LSF; LSF1D; SEF; TFCP2C |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_001166924</u> |
| Locus ID: | 7024 |
| Cytogenetics: | 12q13.12-q13.13 |
| Protein Families: | Transcription Factors |



[View online »](#)