

## Product datasheet for **TS430243P5**

### HABP2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HABP2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230243
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HABP2 Mouse Monoclonal Antibody [Clone ID: OTI8D6] (TA809299)
ACCN:	<u><a href="#">NM_001177660</a></u> , <u><a href="#">NP_001171131</a></u>
Synonyms:	FSAP; HABP; HGFAL; NMTC5; PHBP
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001171131</a></u>
Locus ID:	3026
Cytogenetics:	10q25.3
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)