

Product datasheet for **TS430501P5**

DPP10 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DPP10, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230501
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DPP10 Mouse Monoclonal Antibody [Clone ID: OT11A5] (TA503585)
ACCN:	<u>NM_001178034</u> , <u>NP_001171505</u>
Synonyms:	DPL2; DPPY; DPRP-3; DPRP3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001171505</u>
Locus ID:	57628
Cytogenetics:	2q14.1
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)