

Product datasheet for **TS430884P5**

Dynamin 2 (DNM2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DNM2, transcript variant 5, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230884
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001190716 , NP_001177645
Synonyms:	CMT2M; CMTD11; CMTDIB; DI-CMTB; DYN2; DYNII; LCCS5
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001177645
Locus ID:	1785
Cytogenetics:	19p13.2
Protein Families:	Transcription Factors
Protein Pathways:	Endocytosis, Fc gamma R-mediated phagocytosis



[View online »](#)